

分液漏斗振盪機 Funnel Shaker



VD-12



V-U




VD-12-2S


配件(選購品) / Accessories : (For optional purchases)



離心管固定盤 TH-A
Centrifuge tube holder



試管架固定箱 TH-B
Test tube rack holder



三角燒杯固定箱 FH-B01
Triangular flask holder



1L分液漏斗支持夾
Funnel holder
FH-A



2L分液漏斗支持夾
Funnel holder
FH-B

特點 :

- 本機是將多種試料在同一條件下振盪攪拌均勻，適用於有機與無機化學、醫學、藥品等之基礎研究，或者臨床之溶媒抽出培養產業、排水的水質檢查等各種試驗。
- 指針式振盪次數指示錶，並有定時器。

Feature :

- This shaker is efficient in shaking reagents of a variety of materials under same control parameters. It's suitable for organic/inorganic chemical, medical, pharmaceutical research, extraction and development of media, and waste water analysis.
- This model is with a tachometer and a timer.

Model	100ml	200ml	300ml	500ml	1L	2L	oscillations (epm)	Stoke(mm)	Size (W x L x Hmm)	KG
VS-6	6pcs	6pcs	4pcs	4pcs	3pcs	--	0 ~ 350	40	430 x 430 x 580	59
VD-12	12pcs	12pcs	8pcs	8pcs	6pcs	--	0 ~ 350	40	485 x 430 x 580	65
VD-12-2S	12pcs	12pcs	8pcs	8pcs	6pcs	4pcs	0 ~ 300	40	430 x 710 x 730	81
V-U	6pcs	6pcs	4pcs	4pcs	3pcs	--	0 ~ 350	40	455 x 430 x 580	60

Power : AC 110V/220V, 50/60Hz