

GT SONIC D-series ultrasonic cleaners with degas and double powers



GT SONIC D-series ultrasonic cleaners with digital timer, heater, degas, double powers and memory features. Normal/soft powers control expands the cleaning applications for normal and delicate parts. Degas feature for liquid degassing or beyond cleaning applications.

Accessories:



Basket



Lid



Cord



Features

1. Normal and soft power control to resolve the cleaning blind spot
2. Degas function for improve the cleaning effect.
3. Time Setting: 1-99 minutes adjustable.
4. Heating: 0-80°C temperature adjustable
5. User-friendly and clear panel
6. Memory function
7. LED display for temperature and timer
8. High-performance transducer
9. Patent products

Specifications



Model	Tank size	Unit size	Volume	Ultrasonic power	Ultrasonic frequency	Heating power	Timer	Heating temperature
	L×W×H(mm)	L×W×H(mm)	(L)	Normal (W)	(kHz)	(W)	(MIN)	(°C)
GT SONIC-D2	150×140×100	190×170×220	2	50	40	100	1-99	0-80
GT SONIC-D3	240×140×100	270×170×240	3	100		100	1-99	0-80
GT SONIC-D6	300×155×150	330×180×310	6	150		300	1-99	0-80
GT SONIC-D9	300×240×150	330×270×310	9	200		300	1-99	0-80
GT SONIC-D13	330×300×150	360×330×310	13	300		400	1-99	0-80
GT SONIC-D20	500×300×150	550×330×310	20	400		500	1-99	0-80
GT SONIC-D27	500×300×200	550×330×360	27	500		500	1-99	0-80