

無刷高速型離心機 Brushless High Speed Type



CB-10001

特點：

- 採無碳刷馬達，無污染、防干擾。
- LCD數位面板，具轉速、離心力切換設定顯示功能。
- 10組記憶系統，可預先設定常用參數。
- 不平衡斷電、過載、故障指示安全防護裝置。

Feature：

- Brushless motor. No pollution and no interference
- LCD digital panel with RPM & RCF switch, setting & display functions
- 10 memory programs for pre-setting parameters
- With auto power-off safety devices included rotor imbalance, overload, and error indication

Model	CB-10001
轉速 Speed (rpm)	1,000~10,000
最大離心力 Max. RCF (xg)	10,630
最大容量 Max. Capacity (ml)	400
定時器 Timer (min:sec)	99:59
重量 Weight (kg)	20
尺寸 Size (W x L x H mm)	360 x 480 x 340
電源 Power	AC 110/220V, 50/60Hz

◎ 離心試管另購。 Tubes not included.

■ Rotor sets for CB-10001



RACH-1508



RA-1508-1



RA-1512H-1



RA-5004



RA-5006



RA-5008

Rotor	Capacity (ml x pcs)	Tube Size (Ø x L) mm	Max. Speed (rpm)	Max. RCF (xg)
RACH-1508	15ml x 8	15~17.5 x 92~122	10,000	10,630
RA-1508-1	15ml x 8	15~17.3 x 92~122	8,000	6,660
RA-1512H-1	15ml x 12	15~17.3 x 92~122	6,500	5,010
RA-5004	50ml x 4	27.5~29.5 x 84~117	7,500	6,040
RA-5006	50ml x 6	27.5~29.5 x 84~117	7,000	5,260
RA-5008	50ml x 8	27.5~29.5 x 84~117	6,500	5,010