

WATER BATH

The standard of the water bath

Essential equipment in the laboratory




Perforated heater cover
for safely placing samples
with effective thermal circulation.




WD-11N


WD-23N



Excellent corrosion and chemical resistance stainless steel heater made of **STS #316**



Temp limit device for conducting safe experiments with temperature-sensitive samples



Model : WD-11N

NEW

Water Bath




Features

- Equipped overheat protection sensor.
- Time set function and constant temperature.
- Uniquely designed smooth exterior and mono piece bath.
- Equipped internal bottom for the safety of the user and protection of the unit.

Specifications

Made in Korea

Spec. / Model		WD-11N	WD-23N
Type		Digital	
Controller		PID Controller	
Range		Ambient 5~99 °C	
Temp	Uniformity	±0.5°C	
Size	Internal (W x D x H, mm)	300x240x200	515x330x200
	External (W x D x H, mm)	360x370x280	580x410x280
Safety		Over Temp. Cut-Off	
Material	Internal	Seamless Stainless steel	
	External	Steel, Powder coating	
Electric Requirements		220VAC 50/60Hz	
Heater		1000W	2400W



Model : WD-11N

(mm)								
	Model	W	D	H	W1 (W1')	D1 (D1')	H1 (H1')	g
Dimension	WD-11	360	375	285	300 (285)	240 (225)	200 (165)	35
	WD-23	580	415	281	500 (465)	300 (265)	200 (165)	35

Internal size