

Multi-Position Magnetic Stirrer



ST-500



ST-900



ST-1200

Stainless Steel #304
Speed Adjustable
Separate Control

Model	Position	Dimension (cm)	Speed	Capacity
ST-500LED	5	33 × 36 × 9	MAX 1200 rpm	1000 ml
ST-900LED	9	40 × 42 × 9	MAX 1200 rpm	1000 ml
ST-1200LED	12	44 × 35 × 9	MAX 1200 rpm	1000 ml
ST-1500LED	15	55 × 35 × 9	MAX 1200 rpm	1000 ml