

Portable Pressure Steam Sterilizer / Autoclave



The product utilizes saturated steam to rapidly and efficiently sterilize articles. It can be used in hospitals, public health centers, health stations, clinics in factories and mines, scientific research institutes, etc, to sterilize medical equipment, surgical dressings, glass utensils, solutions, substrata. It can also be used as cooker in plateau areas and be used to make high quality drinking water in industrial and mining enterprises.

Features:

AE-YX18HD

- ◆ Fully stainless steel structure.
- ◆ Hand wheel type of quick-open door structure.
- ◆ Digital display of working status, touch type key.
- ◆ Automatically shut off with beep reminding after sterilization.
- ◆ Electric heated.
- ◆ Easy to operate, safe and reliable.

AE-YX18LD

- ◆ Fully stainless steel structure.
- ◆ Double scale indication pressure gauge.
- ◆ Electric heated.
- ◆ Easy to operate, safe and reliable.

AE-YX18LM

- ◆ Fully stainless steel structure.
- ◆ Electric or LPG heated.
- ◆ Working pressure: 0.14~0.16MPa.
- ◆ Working temperature: 126°C
- ◆ Double scale indication pressure gauge.
- ◆ Easy to operate, safe and reliable.

Technical Specifications:

Model	YX-18HDD	YX-24HDD	YX-18LD	YX-24LD	YX-18LM	YX-24LM
Sterilizing volume	18L (280*260mm)	24L (280*390mm)	18L (280*260mm)	24L (280*390mm)	18L (280*260mm)	24L (280*390mm)
Working pressure	0.14-0.16MPa					
Working	126℃					
Temperature adjust range	105-126℃		126℃			
Timer range	0-99 min		---			
Max. safety pressure	0.165Mpa					
Heat average	≤±1℃					
Power	AC220V.50Hz / 2KW					
Dimension	410×410×630mm	410×410×750mm	410×410×520mm	410×410×630mm	410×410×430mm	410×410×550mm
G.W / N. W	14/12 KG	15/14 KG	18/16KG	19/17 KG	16/14 KG	17/15 KG