

Overhead Stirrer

◎ Large capacity mechanical mixing equipment.

Characteristics:

- Stainless steel head restraints, more durable.
- Digital display model, stirring has three modes: positive, reverse, positive and negative alternating.
- Digital type: brushless motor.
- SH-II-8C is new big power type, high torque, can be used for mixing mud.



Main Technical Parameters:

Patent No. : ZL 2021 3 0287232.0

Model		SH-II-6	SH-II-6C	SH-II-7C	SH-II-8C
Function	Max. Volume (L)	20	20	40	40
	Digital	—	LCD	LCD	LCD
	Speed (r.p.m)	100~2500	100~2500	100~2500	100~2500
	Working Time	Continuous			
	Minimum Resolution (r.p.m)	100	100	100	100
Structure	Exterior	Plastic shell			
	Driving Mode	Motor			
	Motor Power (kW)	0.06	0.06	0.08	0.12
Controller	Continuous Operation	Yes			
	Setting Mode	Scale	Digital	Digital	Digital
	Speed Display Mode	Scale	Digital	Digital	Digital
	Timer	No	Yes	Yes	Yes
	Max. Viscosity(Mpas)	10000	10000	15000	50000
Specification	Exterior Size (W*L*H)(mm)	115 × 95 × 240			
	Packing Size (W*L*H)(mm)	400 × 250 × 210			
	Plate Size (mm)	200 × 300			
	Plate Load	50KG			
	Power Supply (50/60Hz) Current Rating	100~240V			
	NW/GW (kg)	7/8	7/8	7.5/8.5	8/9