

SP-200T Multi-Position Magnetic Hot Plate Stirrer

Magnetic Series

Introduction

SP-200T Multi-position Magnetic Hot Plate Stirrer is designed for the high-efficiency experiments and has 8 stirring bar configurations. It can use laboratory space efficiently. DC brushless motor makes the model, low speed but stable, high speed but powerful with low noise. Stainless steel is coated by silicone resin film disk. Non-slip and anti-corrosion.

Features

1. Digital display stirring speed at real-time, stirring bars operate synchronously and provide uniform agitation reaction conditions.
2. Efficient multi-position magnetic stirrers with 8 stirring positions to use lab space efficiently.
3. DC brushless motor. Low speed but stable, high speed but powerful. Low noise and free maintenance.
4. Stainless steel is coated by silicone resin film pad. Non-slip and anti-corrosion.
5. Applies special heating technology, Surface temperature up to 120 °C.

Parameters

Model	SP-200T(Heating)
Stirring Point Quantity	8
Size of Stirring Bars	Φ8×46mm
Max. Stir Capacity	400ml×8
Digital Display	Temp. & speed
Speed Range	200~1200rpm
Temp. Control Range (surface)	R.T.+5°C~120°C
Platform Material	stainless steel with silicon pad
Motor	DC brushless motor
Voltage	AC 220V/ AC 110V, 50/60Hz
Power	400W
Fuse	250V,5A/3A, Φ5×20
Dimension	W.205×D.480×H.50mm
Net Weight	5.3kgs

