

Water Activity Meter Model WA-60A

It is often used to measure the dryness enduring ability of the microorganisms.

By measuring the water activity of food and selecting reasonable packing and storage methods, the use of preservatives can be reduced and the shelf life of food, grain and fruits and vegetables can be judged.



Sensor	Non-conductive humidity sensor (Stainless steel)
Range	0 – 1.0 Aw
Resolution	±0.01 Aw
Accuracy	±0.02 Aw
Sampling Time	5 minutes
Display	Backlight LCD (Aw & Temperature)
Operating Conditions	0–50°C, <95%RH
Power Supply	2 × 1.5V AAA batteries