

# CAPP

D E N M A R K

## Capp *Bravo*

SINGLE CHANNEL PIPETTES

Manual pipetting has never been easier! Capp Denmark has combined lightweight ergonomic design with the latest innovations in manual liquid handling to create Capp *Bravo*, one of the most user friendly manual pipettes on the market today.

Capp *Bravo* is constructed of the highest quality autoclavable materials and features feather light plunger and ejection forces. There is a built-in oversized indicator for easy volume monitoring and precise adjustment. With Capp *Bravo* all of these features are packaged in a modern, streamlined, world-renowned Danish design. Experience Capp's symphony of precise measurements and see for yourself why it is called *Bravo*!

### Safety:

Fully Autoclavable  
Positive indicator lock

### Ergonomics:

Ergonomically constructed hook for comfortable handling  
Smooth tip mounting and ejection

### Controls:

Oversized indicator  
Reduced force feather touch plunger and ejector

### Unique Features:

Metal shaft with double o-rings for ultra micro sizes  
Ultra lightweight construction





 CAPP APS  
 LANDBRUGSVEJ 10  
 DK-5260 ODENSE  
 DENMARK

TEL. (+45) 6613 6140  
 FAX. (+45) 6613 2770

MAIL: INFO@CAPP.DK  
 WEB: WWW.CAPP.DK



## Ordering Info

### Capp *Bravo* variable volume

| Cat. No. | Description                                                        | Inaccuracy % | Imprecision % | Color                                                                                 |
|----------|--------------------------------------------------------------------|--------------|---------------|---------------------------------------------------------------------------------------|
| B02-1    | Capp pipette, variable vol. 0.1-2 ul, Bravo                        | 12.00 / 2.00 | 6.00 / 1.00   |  |
| B10-1    | Capp pipette, variable vol. 0.5-10 ul, Bravo                       | 2.50 / 1.00  | 1.80 / 0.50   |  |
| B20-1U   | Capp pipette, variable vol. 2-20 ul, for use with 10ul tips, Bravo | 5.00 / 1.00  | 1.50 / 0.30   |  |
| B20-1    | Capp pipette, variable vol. 2-20 ul, Bravo                         | 5.00 / 1.00  | 1.50 / 0.30   |  |
| B50-1    | Capp pipette, variable vol. 5-50 ul, Bravo                         | 2.00 / 1.00  | 1.20 / 0.30   |  |
| B100-1   | Capp pipette, variable vol. 10-100 ul, Bravo                       | 2.50 / 0.70  | 1.00 / 0.20   |  |
| B200-1   | Capp pipette, variable vol. 20-200 ul, Bravo                       | 2.50 / 0.60  | 0.70 / 0.20   |  |
| B1000-1  | Capp pipette, variable vol. 100-1000 ul, Bravo                     | 2.00 / 0.80  | 0.60 / 0.20   |  |
| B5000-1  | Capp pipette, variable vol. 1-5 ml, Bravo                          | 1.20 / 0.60  | 0.30 / 0.15   |  |
| B10000-1 | Capp pipette, variable vol. 1-10 ml, Bravo                         | 3.00 / 0.60  | 0.60 / 0.15   |  |

### Capp *Bravo* fixed volume (customized volumes available on request)

| Cat. No.  | Description                             | Inaccuracy % | Imprecision % |
|-----------|-----------------------------------------|--------------|---------------|
| B01-1F    | Capp pipette, fixed vol. 1 ul, Bravo    | 2.50         | 1.80          |
| B02-1F    | Capp pipette, fixed vol. 2 ul, Bravo    | 2.00         | 1.00          |
| B05-1F    | Capp pipette, fixed vol. 5 ul, Bravo    | 1.20         | 0.80          |
| B10-1F    | Capp pipette, fixed vol. 10 ul, Bravo   | 1.00         | 0.50          |
| B20-1 F   | Capp pipette, fixed vol. 20 ul, Bravo   | 1.00         | 0.30          |
| B50-1F    | Capp pipette, fixed vol. 50 ul, Bravo   | 1.00         | 0.30          |
| B100-1F   | Capp pipette, fixed vol. 100 ul, Bravo  | 0.70         | 0.20          |
| B200-1F   | Capp pipette, fixed vol. 200 ul, Bravo  | 0.60         | 0.20          |
| B500-1F   | Capp pipette, fixed vol. 500 ul, Bravo  | 0.80         | 0.15          |
| B1000-1F  | Capp pipette, fixed vol. 1000 ul, Bravo | 0.80         | 0.20          |
| B5000-1F  | Capp pipette, fixed vol. 5 ml, Bravo    | 0.60         | 0.15          |
| B10000-1F | Capp pipette, fixed vol. 10 ml, Bravo   | 0.60         | 0.15          |

### Capp *Bravo* Starter Packs

| Cat. No. | Description                                                            |
|----------|------------------------------------------------------------------------|
| SK-03-BB | Three pipettes free of choice, C-05 stand, 3 racks of tips and 1 timer |

### Pipette stands

| Cat. No. | Description                        |
|----------|------------------------------------|
| C-05     | Pipette stand for up to 5 pipettes |