

Centrifuge

Lugar

model LC-04B



Introduction :

The frame of this instrument is made of high strength plastic. Its model is beautiful, and it has the advantages of small cubage, low weight, large capacity, low noise, high efficiency and so on. It can be used in hospitals and biochemical labs for qualitative analysis of serum, urea and plasma. Rotor imbalance protection device, speed choice in advance & acceleration accurate gradually, door switch

Technical Parameters :

- Max. Speed : 4,000 r/min
- Max. RCF : 2325xg
- Capacity : 12 x 20 ml
- Time range : 0 - 30 min
- Power source : 220V 50Hz /110V 60Hz 135W
- Dimension : 280*310*265(L*W*H)
- Packing Size : 36*34*35(cm)
- Net weight : 10Kg
- Products Made in China
- Warranty 1 year