

# Centrifuge

**Lugar**

**model LC-04C**



## Introduction :

- ◆ Strong Plastic case with metal lid
- ◆ High Brushless motor, with mechanical door lock
- ◆ Over speed and imbalance protection
- ◆ With safety switch, stop working when open the door
- ◆ Stepless speed adjustment , LED display speed and timer
- ◆ Error code display, speed up and down adjustment

## Technical Parameters :

- Max. Speed : 4,000 RPM
- Max. RCF : 2325xg
- Capacity : 12 x 20 ml
- Time range : 0 -99 min
- Power source : 220V 50Hz
- Dimension : 280\*310\*265(L\*W\*H)
- Weight : 9Kg
- Products Made in China
- Warranty 1 year