

Refrigerated Heating Circulator

HS-3005HC

As a source of cooling water to the Rotary Evaporator and as a heat-emitting equipment it supplies sufficient quantity of cooling water as highly efficient chemical pumps are installed and does not contaminate other instruments as it does not get corroded.



MODEL	HS-3005HC
Capacity	7 Liter
Bath Dimension (W × D × H)	220 × 270 × 145 mm
Dimension (W × D × H)	430 × 350 × 540 mm
Materials	Stainless Steel
Heater	800 W
Temp. Range	-20 °C to + 80 °C
Temp. Sensor	PT 100
Circulation Pump	15 W, 15 Liter/min
Controller	FND
Refrigerator	500 W
Refrigerant	R404
Gross Weight	31.3 Kg
Power Supply	AC 230 V, 50/60 Hz