

Rotary Evaporator

HS-5SP (5 Liter)

Distinguished Features


- Evaporating flask : 5,000 mL
- Receiving flask : 3,000mL
- Digital speed control
- Digital R.P.M Display/Speed control range : 20~280 R.P.M
- Vertical type condenser
- Digital interior temperature measure and display
(It catches the temperature changes caused by the interior pressure drops)
- Motorized UP-DOWN control switch
- Easy to utilize standard distillation solution
- Auto timer function
- Autonomous condenser lifting function after a complete shut down of the machine

Included Water Bath HS-3001

(PTFE Coating)

If the amount of water goes down under the certain threshold in the heating bath, the system automatically turns off its power for the safety purpose

- Included vacuum controller
- *Vacuum range : 0~760 mmHg



Registration No.
20-0461039/30-0668181

HS-5SP (5 Liter)

Heater	1,800 W
Bath Material	SUS 304
Weight	3.5 Kg
Bath Dimension (Ø x H)	300 x 150 mm (Capacity : 8 Liter)
Required Power Source	AC 220 V, 3.5 A, 50/60 Hz
Temp. Control Range	Amb. ~180° C
Temp. Setting/Display	Digital Control (Controlled by the main controller)

*Double jacketed bath body (Safe to touch)

MODEL	HS-5SP (5 Liter)
Dimension (W x D x H)	400 x 400 x 800 mm
Weight	26 Kg
Power Source	AC 220 V, 2.5 A, 50/60 Hz