

D1012UA

High Speed Palm Micro Centrifuge

Applications

This centrifuges is ideal for quick spins,microfiltration of samples,cell separations and many other routine laboratory procedures.

Features

- Max.Speed 10000rpm (7490xg) ;
- OLED display, may switch speed and time;
- DC motor drive, smooth and quiet operation;
- Safe,endurable and robust rotor with clamp locking design, easy to operate;
- Wide voltage technology: the speed is not affected by voltage fluctuations, high speed accuracy;
- Compatible with 12×0.5/1.5/2mL tubes+4×0.2mL PCR 8-tube strip, 44×0.2mL PCR tubes;
- Screen buttons can adjust the start/stop of the motor, supports both soft braking and instant stop modes;
- Simple appearance and structure, only weight 1.5Kg;



Specifications	D1012UA
Speed Range	2000-10000rpm,increment:1000rpm
Max.RCF	6708×g-1.5/2mL; 6596-7490×g-PCR8 strips;
Rotor Capacity	44×0.2mL PCR tubes 12×0.5/1.5/2mL tubes+ 4×0.2mL PCR 8-tube strip
Run Time	10s-20min
Display	OLED
Acceleration/Deceleration	1 ↑ / 2 ↓
Driving Motor	DC motor
Power	AC100-240V/50Hz/60Hz 40W
Noise Level	≤58dB
Dimension (D × W × H)	210×180×130mm
Weight	1.5kg

