



BJ-2 DISINTEGRATION TESTER

Features

There are 12 tubes in two basket assembly.

Each test station is driven independently. Auto stop at timing point. LCD display.

The bath liquid can achieve an even temperature.

A MPU is used to automatically control the temperature, time and frequency of nacelle back and forth.

Timing at will. The preset and real time data can be displayed alternately in time.

Automation: auto-test, auto-diagnose, auto-alarm

Specifications

Basket assembly: 2

Nacelle back and forth frequency: 30~32 /min

Nacelle back and forth range: 55±1 mm

Time range: 1~900 min

Temperature range: RT~45°C

Temperature accuracy: ±0.5°C

Power: 220V/ 50Hz/ 600W or 110V/ 60Hz/ 600W

Dimension: 400×300×420 mm