




Features

- ▶ High speed mixer for the research & experiment purpose of emulsification dispersion, homogeneity and dissolution
- ▶ Suitable in various fields such as medical, science, cosmetics and foods
- ▶ High speed disperser for using the uniform dispersion ideally and easy to exchange the stirrer blade
- ▶ Speed of 1000~12000 rpm
- ▶ Suitable for the sample of low viscosity
- ▶ Constant speed by feedback control
- ▶ Digital display

*Fixing clamp (option)



Specifications

Model	HM1200D	HD1200D
Motor(output)	142W	
Working Capacity	500~2000ml	500~5000ml
Speed Range	1000~11000rpm	1000~12000rpm
Speed Controller	Feedback control	
Speed Display	Digital	
Dimensions (W×D×H)	70×180×170mm	
Weight	3.2Kg	2.8Kg
Electrical Supply	AC 220V 50/60Hz	
Price		
Price (Set)		

Note : Set ⇒ Mixer or Disper + Stand(PS80) + Clamp(CH320)