

TECHNICAL SPECIFICATION

MODEL	LC-MSB-S	LC-MSB-H	LC-MSB-D	LC-MSB-HD
Voltage (V)	100-240	110/220	100-240	110/220
Powe (W)	10	450	24	550
Display Method	/	/	LCD	LCD
Working Plate Material	Tempered Glass	Cast aluminum +Ceramic spraying	Tempered Glass	Cast aluminum +Ceramic spraying
Speed Range (rpm)	200-1600		100-1600	
Max Stirring Volume(L)	5L			
Timing	/	/	1-99h59min	1-99h59min
Heating Function	NO	YES	NO	YES
Temprature Control Range	/	RT-300°C	/	RT-300°C
Temperature Control Accuracy	/	/	/	Working plate±3 Liquor±0.5
Working Plate Size(MM)	150*150			
Motor Type	DC Brushless Motor			
PT1000	/	/	/	Standard
Program Function	/	/	YES	YES
Use Environment	5-40°C 80%RH			
Dustproof and waterproof rating	IP42			



LC-MSB-S



LC-MSB-H



LC-MSB-D



LC-MSB-HD