

TECHNICAL SPECIFICATION

MODEL	LC-OES-60	LC-OES-120	LC-OES-200	LC-OES-60SH	LC-OES-120SH	LC-OES-200SH
Voltage	AC100-240					
Frequency (Hz)	50-60					
Power (W)	72	120	200	72	120	200
Display Mode	LCD Display					
Max Stirring H ₂ O Volume (L)	20	40	60	20	40	60
Suggest Stirring H ₂ O Volume (L)	10	15	20	10	15	20
Max Stirring Viscosity (mPa·s)	10000	20000	30000	10000	20000	30000
Suggest Stirring Viscosity (mPa·s)	3000	7000	10000	3000	7000	10000
Speed Range (rpm)	100-2000					
Timing Range	/			1-99h59min		
Max Torque (N·cm)	10	20	30	10	20	30
Motor Type	DC Permanent Magnet Motor					



LC-OES-60/120/200



LC-OES-60SH/120SH/200SH

MODEL	Standard Mixing Propeller	Specifications
LC-OES-60	Two blade mixing propeller*1	250*6*45mm
LC-OES-120		300*8*65mm
LC-OES-200		300*8*65mm
LC-OES-60SH	Two blade mixing propeller*1 Three blade mixing propeller*1	250*6*45mm / 250*6*30mm
LC-OES-120SH		300*8*65mm / 250*8*50mm
LC-OES-200SH		300*8*65mm / 300*8*65mm